

QUANTIFIT®

INTRODUCING...FITTRACK GOLD™

FitTrack Gold offers a completely new experience for Quantifit users. The software has an updated interface that's easier to use, including many tools that will revolutionize the way data is managed. With a SQL Express database, users will experience more flexibility in data storage as well as reporting functions. Here are some of the exciting features in the new FitTrack Gold:

The screenshot shows the FitTrack Gold 2.0 dashboard with the following sections and callouts:

- Navigation:** A sidebar menu on the left includes Home, Personnel, Company, Reports, Tools, Help, Admin, and Exit. A callout points to this menu with the text "Easy Navigation".
- Due for Test:** A table listing personnel with columns for Personnel ID, Name, and Job Code. A callout points to the first row with the text "See who's due for a test".
- Links:** A section containing web links for OHD, Service Return Form, and OSHA 1910.134 App. A callout points to this section with the text "Convenient web links".
- Mask Inventory:** A table showing mask models, sizes, and counts. A callout points to this table with the text "Mask inventory currently in use".
- Calibration Dates/Due:** A table listing Quantifit IDs and their NIST Calibration Due Dates. A callout points to this table with the text "Keep track of instrument calibration dates".
- Test Statistics:** A table showing test results by week, month, and year. A callout points to the first row with the text "Test tracking at a glance".
- Quantifit Detail:** A section at the bottom right with a message "An attached Quantifit was not found." and a "Connect" button.

DASHBOARD

A new dashboard home screen displays useful information at a glance. This includes a list of people who are due to be tested, test statistics to track the number of tests performed, a mask inventory for the current year of testing, calibration due dates for your instruments, and web links that can be customized to your specific needs.

SECURITY

Security is comprehensive and flexible. User names and passwords can be assigned, and the administrator can control any task within the software.

UPGRADING FIRMWARE

Instrument updates can be performed directly through the software.